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
Product introduction

Health starts with teeth. Thank you for choosing our LED household intelligent ultrasonic household tooth cleaner. The LED lamp household intelligent ultrasonic household tooth cleaner is specially designed for family users. Through innovative design, the complex structure of the original tooth cleaner is changed into an intelligent chip integrated circuit, and the advanced intelligent frequency stabilization design concept is adopted. It has stable performance and strong applicability, and can cooperate with a variety of functions. The whole process of tooth cleaning is completed under the control of microprocessor, which ensures the stable output of power and can effectively remove dental calculus, smoke stains, plaque and tea stains. The operation is simple, convenient, fast and efficient.

How to use it safely

- If you use any improper method that is not described in this instruction manual or is not used correctly, please stop using it. In case of failure, please stop using.
- This product is not a toy. Please don't let your child use it directly.
- It is forbidden for patients with pacemaker.
- Patients with heart disease, pregnant women and children should be cautious.
- Please do not put the product in direct sunlight, it should be stored in 10° to 40°C indoor environment.
- Please don't disassemble, transform or bend this product by yourself.
- In case of product quality problems, please contact the customer service of the store where you purchased this product.

Product features

1. This product is designed for family, simple, fast and efficient.
2. UHF vibration, gentle, comfortable, enhanced three mode adjustable.
3. 10 minutes automatic stop work protection design, stop work after the LED light off.
4. Integrated design, head LED lighting, clear observation of oral conditions.
5. The ultrasonic transducer is made of aviation titanium alloy.

Product principle

The high-frequency electric pulse wave is generated by the electronic oscillation circuit, amplified to a certain intensity, and transmitted to the transducer, which converts the electric energy into mechanical energy, and stimulates the working tip to produce the same frequency vibration, so as to achieve the purpose of crushing the dental calculus and sweeping away the stubborn dental stains and loose tartar. The ion is produced by the ion generator, and the ion is produced by holding the ion connecting ring. The negative ions are transmitted to the teeth, and the tooth stains are loose and separated from the teeth by using the principle of negative ion exclusion.

Matters needing attention

1. When the ultrasonic working tip touches the teeth, it makes a beeping sound, which is the normal working state.
2. If you have periodontal disease, consult your dentist before use.
3. If it is on the boundary of the gingiva, gingival bleeding may occur even in the case of healthy gingiva.
4. If gum bleeding persists, which may be due to gum defects such as periodontal disease, discontinue use and consult your dentist.
5. In the process of tooth washing, do not stay on the surface of teeth for a long time with the working tip. To remove the calculus and stains of stubborn teeth, please use the working tip to remove them.
6. When working without water, do not touch the gums for a long time. Long time contact with the gums may bum the gums.
7. After the machine is started, there is no peak sound at the working point. If there is a peak sound, it is a normal phenomenon, you can re install the working tip or tighten the working tip.
8. Even if it is not used directly with the front end of the head, a weak current will be felt due to the ion output when treating the area near the parts with tooth decay, gold crown or silver crown. In this case, please do not touch the ion connection ring, only use the ultrasonic output for treatment, and the ion output will not be transmitted to the teeth.


Cleaning and maintenance

1. Do not put or drop the product into or into other liquids.
2. The equipment should be handled with care and stored in a dry place.
3. Do not mix with toxic, corrosive, flammable and explosive materials.

Matters needing attention

1. Do not soak in water or liquid.
2. Keep the machine clean before and after use.
3. The working tip must be tightened when using.
4. Power on for the first time to use, first run for 2 ~ 3 seconds, stop for 1 second, such operation several times, easy to work head in the best state.
5. Do not bend or grind the work tip.
6. Do not knock or scrape the handle.
7. Please dry the machine after use and keep it in a dry and cool environment.
8. Do not pull the power cord during the use of the equipment, so as not to cause the handle to separate from the machine.


Popis produktu


How to operate

1. Make sure the cleaner is fully charged before use.
2. Press the on key to start the machine, turn on the light and enter the soft gear.
At the same time, turn on the negative ion and ultrasonic mode. Press the on key again, the negative ion remains unchanged, and the ultrasonic enters the comfort gear. And so on, there are three gears in total. It is suggested to start from the low gear and gradually increase the gear after adaptation.
 1. Soft
 2. comfort
 3. strong
3. Select a suitable working point and install it clockwise on the main engine with a small wrench. DO not use strong rotation, so as not to damage the interface. First, open the appropriate gear and vibrate the tooth stone with the metal head. Then, when operating on the tooth stone, scrape it with a little strength and try more times. Because the tooth stone is stubborn, a little more patience will make it easier to master the method.
4. Press C on and off, the mode light shines, and the LED light at the front end of the handle lights up, indicating that the machine is running.
5. When the machine works normally, the vibration frequency is very fast. When cleaning teeth, it only needs light contact and repeated movement at a certain speed to clear teeth stones.
Caution: DO not apply partial force or stay too long when cleaning teeth!
6. When off, please press the on key, and unplug the power plug after the LED light is off.
7. The machine is a precision equipment. Please avoid damage to internal and external components caused by fling, touching and dropping.
8. Brush head is to clean the teeth after calculus, can not be used as an electric toothbrush.
9. Low voltage detection, when the battery voltage is lower than 3V, flash light alarm.
10. After use, long press the open key, the power off, the equipment preset 10 minutes automatic shutdown function, to prevent overheating protection.

Introduction to work tips

Tip 1 is thin, used to remove calculus and tooth smoke stains on gingiva.


Tip 2 is wider, used to remove large stones and smoke stains on the neck and lingual side of the supragingival teeth.

- Please steam, soak in alcohol and disinfect the working tip by high temperature and high pressure before use.
- In order to prevent cross infection, do not share the work tip with many people.
- The machine frequency is different, please choose the right working point.
- Please use the work tip correctly as shown above.


How to use work tips

The work tip gently touches the teeth vertically and moves to the top and bottom as well as left and right. Please adjust the use time of each part according to the condition of teeth.


⚠ Please do not use the surface of the work tip without tooth dirt.


Tooth surface dirt


Boundary between gum and tooth


Interdental gap



Dirt from uneven teeth

Dirt on the back of teeth

Please use mouth glasses when treating the back of teeth


How to change the work tip


The orientation of the work tip may be different but as long as it is firmly mounted on it, it can be used in any direction.


If the connection is loose, the output cannot be transmitted correctly, so make sure the work tip is tightly connected and use the wrench kit for the work tip installation. Please do not use other tools to disassemble the sample plate, which may cause failure or damage.


Confirm position A

Insert the wrench in position A

Rotate the head counter clockwise with a wrench to separate the working tip and complete the disassembly.


Preventive measure

Avoid over tight assembly, resulting in damage to the interface.

How to maintain

Please do not use diluents, benzene and solvents, chemicals, scrubbing powder.

The main body can be directly washed with water, and the dirt on the main body can be wiped with a rag.